EXHIBIT P THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE

Agency		GOODYEAR Date 1/23/2004 Time 1204
		erial # 66-06 4416 Location GP6
QA Specialist Wile		
(Print Name)		
(5)	1.	Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0/66 AC
$(\boldsymbol{\nu})$	2.	Turn on the simulator and allow the temperature to reach 34°C±.2°C.
(~)	3.	Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
(_ ~)	4.	Attach simulator to the simulator entrance port on the Intoxilyzer.
(2)	5.	Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
(,u)	6.	Push Start Test button or press enter on keyboard.
AL STRETE	7.	Insert card in response to display.
(2)	8.	Air blank completed.
(Ju)	9.	Calibration check completed. Test results 0. 295 AC
(10.	Air blank completed.
(,,)	11.	When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
(\sigma)		Return mode selection switch #9 to off position after all calibration checks are complete or type Q and
enter on keyboard.		
SIGNATURE		
*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD		

DHS/DSLS/Form C116 (Rev. 7-93)

GOODYEAR POLICE DEPT.
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
01/23/2004 12:04
COPY NUMBER 1 OF 2

OFFICER BADGE=0000974
N.C.I.C NO.=AZ0071500

TEST BRAC VALUE TIME
AIR BLANK .000 12:03
CAL. CHECK .095 12:04
AIR BLANK .000 12:04

OPERATOR:

ADDITIONAL INFORMATION: