## EXHIBIT P THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C R9-14-404(A)

#### ARIZONA DEPARTMENT OF HEALTH SERVICES

# STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

#### STANDARD CALIBRATION CHECK PROCEDURE -

Agency GODDY EAR P.D. Date 6-13-2004 Time 2/29				
Intoxilyzer Serial # 6.6-00 4416 Location JUSTICE FACILITY				
QA Specialist Ti Vallandiashan 981				
(Print Name)				
1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the				
simulator. Insure that a tight seal has been made. Standard value: 0/00 AC				
2. Turn on the simulator and allow the temperature to reach 34°C±.2°C.				
3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or				
selecting "C" on keyboard menu.				
( ) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.				
5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".				
6. Push Start Test button or press enter on keyboard.				
7. Insert card in response to display.				
8. Air blank completed.				
( ) 9. Calibration check completed. Test results 0.094AC				
( ) 10. Air blank completed.				
( ) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.				
12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and				
enter on keyboard.				
SIGNATURE TO CONTROL TO TO TO THE SIGNATURE				
*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD				
DHS/DSLS/Form C116 (Rev. 7-93)				
DHS/DSLS/Form C116 (Rev. 7-93)				

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 21:30 06/13/2004 COPY NUMBER 2 OF 2 OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 BRAC VALUE TIME TEST .000 21:29 .094 21:29 AIR BLANK CAL, CHECK AIR BLANK .000 21:30 OPERATOR: ADDITIONAL INFORMATION #

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/13/2004 21:30 COPY NUMBER 1 OF 2 OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK .000 21:29 CAL. CHECK .094 21:29 21:30 AIR BLANK .000 951 OFERATOR: ADDITIONAL INFORMATION:

#### EXHIBIT Q THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R9-14-404(A)

#### ARIZONA DEPARTMENT TOF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

AGENCY (1000YEAR P.D. QA SPECIALIST J. Vallanding 44)
AGENCY (10004EAR P.1). QA SPECIALIST J. VG//4nd.ing 44.  INTOXILYZER SERIAL #66 0044/6 LOCATION JUSTICE FACILITY  DATE 6/3-2004 TIME 2/31
Display reads "READY TO START" or "PUSH BUTTON TO START TEST."
DIAGNOSTIC TESTS  DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is
OPERATIONAL TESTS  Alcohol-free subject Test result 0. OOO AC  Error recognition logic system functioning  Invalid test printed  Proper sample recognition system  Invalid sample printed  Deficient sample printed  Completeness of sample purge with collection tube **  Calibration standard 0. / OO AC Results: 0. O95  Instrument operating properly and accurately. YESNO
COMMENTS: DURING THE ERROR RECOGNITION COGIC SYSTEM (RFI TEST) THE DISPLAY READ "AMRIENT FAILURE" AND SHOULD HAVE READ "INHIBITED -RFI"
SIGNATURE  SIGNATURE  Alcohol concentration standard collection for subsequent analysis.**  *WITH OR WITHOUT VAPOR RECICULATION AND WITH OR WITHOUT KEYBOARD

\*\*THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR SAMPLE CAPTURE

DHS/DSLS/Form C117 (Rev. 7-93)

GOODYEAR PD-JUSTICE FACILITY
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
06/13/2004
COPY NUMBER 1 OF 1

21:31

OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500

DVM TEST  $,2/\rho$ 

OPERATOR:	JUL	7551
ADDITIONAL	INFORMATION:	

GOODYEAR PD-JUSTICE FACILITY
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
06/13/2004
COPY NUMBER 1 OF 1

21:31

OFF	ICER	BADGE=0000981
N.C	.I.C	NO.=AZ0071500

DISPLAY TEST

OPERATOR: _	Such	951
ADDITIONAL	INFORMATION:	

DATE 06/13/2004 TIME 21:33 COPY NUMBER 1 OF 1

OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500

ABCDEFGHIJKLMNOPQRSTUVWXYZ0123456789!@#\$

OFERATOR: 99/
ADDITIONAL INFORMATION:

GOODYEAR PD-JUSTICE FACILITY
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
06/13/2004
COPY NUMBER 1 OF 1

21:33

OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500

TIME, DATE, LOC SET

OPERATOR:	Me	GH 481
ADDITIONAL	INFORMATION	
***************************************		

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/13/2004 21:40 COPY NUMBER 1 OF 2 SUB NAME=TEST, SQAP-AF, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =111 OFFICER =VALLANDINGHAM, J, OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK .000 21:40 SUBJECT TEST .000 21:40 AIR BLANK .000 21:40 OPERATOR: ADDITIONAL INFORMATION:

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/13/2004 21:40 COPY NUMBER 2 OF 2 SUB NAME=TEST, SQAP-AF, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =111 OFFICER =VALLANDINGHAM, J, OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK .000 21:40 .000 SUBJECT TEST 21:40 21:40 AIR BLANK  $_{\bullet}$ OPERATOR: ADDITIONAL INFORMATION:

SN 66-004416 DATE 06/13/2004 COPY NUMBER 1 OF 2

VERSION E 773.23 TIME 21:42

INVALID TEST

SUB NAME=TEST,SQAP-PB,
SUB DOB =11/11/11
SUB SEX =MALE WEIGHT =111
OFFICER =VALLANDINGHAM,J,
OFFICER BADGE=0000981
N.C.I.C NO.=AZ0071500

PRESSED START TEST AT IMPROPER TIME

OPERATOR: _	Silvell 981
ADDITIONAL	INFORMATION:

SN 66-004416 DATE 06/13/2004 COPY NUMBER 2 OF 2

VERSION E 773.23 TIME 21:42

INVALID TEST

SUB NAME=TEST,SQAP-PB,
SUB DOB =11/11/11
SUB SEX =MALE WEIGHT =111
OFFICER =VALLANDINGHAM,J,
OFFICER BADGE=0000981
N.C.I.C NO.=AZ0071500

PRESSED START TEST AT IMPROPER TIME

OPERATOR: _	Sa	Vill 981	
ADDITIONAL	INFORMATION:		-

SN 66-004416 DATE 06/13/2004 COPY NUMBER 1 OF 2

VERSION E 773.23 TIME 21:43

INVALID TEST

SUB NAME=TEST, SQAP-BWT,
SUB DOB =11/11/11
SUB SEX =MALE WEIGHT =111
OFFICER =VALLANDINGHAM, J,
OFFICER BADGE=0000981
N.C.I.C NO.=AZ0071500

SAMPLE INTRODUCED AT IMPROPER TIME.

OPERATOR: _	9.	V S	981
ADDITIONAL	INFORMATION:		
			······································

SN 66-004416 VERSION E 773.23 DATE 06/13/2004 TIME 21:43 COPY NUMBER 2 OF 2

INVALID TEST

SUB NAME=TEST, SQAP-BWT, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =111 OFFICER =VALLANDINGHAM, J, OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500

SAMPLE INTRODUCED AT IMPROPER TIME.

OPERATOR:	Will ass	
ADDITI <b>ON</b> AL	INFORMATION:	
		*************************************

SN 66-004416 VERSION E 773.23 DATE 06/13/2004 TIME 21:45 COPY NUMBER 1 OF 2 INVALID TEST SUB NAME=TEST, SQAP-RFI, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =111 OFFICER =VALLANDINGHAM.J. OFFICER BADGE=0000981 N.C.I.C NO. = AZO071500 SHOULD HAVE SAID CHECK AMBIENT CONDITIONS "INHIBITED -RFI" OPERATOR: ADDITIONAL INFORMATION:

SN 66-004416 VERSION E 773.23 TIME 21:45 DATE 06/13/2004 COPY NUMBER 2 OF 2 INVALID TEST SUB NAME=TEST, SQAP-RFI, SUB DOB =11/11/11 WEIGHT =111 SUB SEX =MALE OFFICER =VALLANDINGHAM, J, OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 SNOULD HAVE SAID CHECK AMBIENT CONDITIONS " INHIBITED - R FI" ADDITIONAL INFORMATION:

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 21:51 06/13/2004 COPY NUMBER 1 OF 2 SUB NAME=TEST, SQAP-MOA, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =111 OFFICER =VALLANDINGHAM, J, OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME .000 21:50 AIR BLANK INVALID SAMPLE .MOA 21:50 AIR BLANK .000 21:51 BEGIN NEW DEPRIVATION PERIOD OPERATOR: ADDITIONAL INFORMATION:

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/13/2004 21:51 COPY NUMBER 2 OF 2 SUB NAME=TEST, SQAP-MOA. SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =111 OFFICER =VALLANDINGHAM, J. OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK .000 21:50 INVALID SAMPLE .MOA 21:50 AIR BLANK .000 21:51 BEGIN NEW DEPRIVATION PERIOD OPERATOR: ADDITIONAL INFORMATIONS

INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/13/2004 21:56 COPY NUMBER 1 OF 2 SUB NAME=TEST, SQAP-DS, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =111 OFFICER =VALLANDINGHAM, J, OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK .000 21:53 \*SUBJECT TEST .012 21:56 AIR BLANK  $"\bigcirc\bigcirc\bigcirc$ 21:56 \* DEFICIENT SAMPLE - VALUE PRINTED WAS HIGHEST OBTAINED. ADDITIONAL INFORMATION:

GOODYEAR PD-JUSTICE FACILITY

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/13/2004 21:56 COPY NUMBER 2 OF 2 SUB NAME=TEST.SQAP-DS, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =111 OFFICER = VALLANDINGHAM, J, OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK 21:53 \*SUBJECT TEST .012 21:56 .000 21:56 AIR BLANK \* DEFICIENT SAMPLE - VALUE PRINTED WAS HIGHEST OBTAINED. OPERATOR: ADDITIONAL INFORMATION:

### EXHIBIT P THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C R9-14-404(A)

### ARIZONA DEPARTMENT OF HEALTH SERVICES

### STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

#### STANDARD CALIBRATION CHECK PROCEDURE -

Agency Gospyear P.D. Date 6/3-2004 Time 2/58 Intoxilyzer Serial # / 66 oc 44/6 Location, Justice Facility  QA Specialist  J. Wallanding Gen. 9ff				
Intoxilyzer Serial # / 66 oc 44/6 Location, JUSTICE FACILLY				
OA Specialist J. Hallandias 46m 9Al				
(Print Name)				
1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the				
simulator. Insure that a tight seal has been made. Standard value: 0. 100 AC				
2. Turn on the simulator and allow the temperature to reach 34°C±.2°C.				
3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.				
4. Attach simulator to the simulator entrance port on the Intoxilyzer.				
( ) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".				
6. Push Start Test button or press enter on keyboard.				
( ) 7. Insert card in response to display.				
( ) 8. Air blank completed.				
9. Calibration check completed. Test results 0.095AC				
( ) 10. Air blank completed.				
( ) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.				
( ) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and				
enter on keyboard.				
cinci on keyoodid.				
SIGNATURE 751				
SIGINITORE				
*WITH OR WITHOUT VAPOR RECERCULATION AND WITH OR WITHOUT KEYBOARD				
DUCTOSI S/Com C116 (Bar. 7.02)				

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/13/2004 21:59 COPY NUMBER 2 OF 2 OFFICER BADGE=0000981 N.C.I.C NO.=AZO071500 TEST BRAC VALUE TIME AIR BLANK .000 21:58 CAL, CHECK 21:58 .095 AIR BLANK  $_{*}\bigcirc\bigcirc\bigcirc\bigcirc$ 21:59 OPERATOR: ADDITIONAL INFORMATION:

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/13/2004 COPY NUMBER 1 OF 2	21:59
OFFICER BADGE=0000981 N.C.I.C NO.=AZ0071500	
TEST BRAC VALUE TIME AIR BLANK .000 21:58 CAL. CHECK .095 21:58 AIR BLANK .000 21:59  OPERATOR:	
ADDITIONAL INFORMATION:	

### EXHIBIT P THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C R9-14-404(A)

#### ARIZONA DEPARTMENT OF HEALTH SERVICES

### STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

STANDARD CALIBRATION CHECK PROCEDURE

Agency_	GOOD GOT	PD	Date	6/14/04	Time //35				
Intoxilyz	Serial # (06 C	04416	Locatio	n <i>GPD</i>					
QA Specialist Will									
(Print Name)									
el.		1 . 1 . 1							
(W)					tor jar and assemble the				
.11.				Standard value: 0. /	<u>U)</u> AC				
• • • • • • • • • • • • • • • • • • • •		Turn on the simulator and allow the temperature to reach 34°C±.2°C.							
(U)	-	Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or							
1,	selecting "C" on keyboard menu.								
(N)		Attach simulator to the simulator entrance port on the Intoxilyzer.							
(W) :	- · · · · · · · · · · · · · · · · · · ·	Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".							
$(\mathcal{W})$	Push Start Test butt	Push Start Test button or press enter on keyboard.							
(N) :	Insert card in respo	ise to display.							
(W) 8	Air blank complete								
(1L) 9	Calibration check c	ompleted. Test resul	Its <u>0.( 0()</u>	AC.					
$(\mathcal{L})$ 1	). Air blank complete	l.	-						
$(\mathcal{U})$ 1	. When display reads	Test Complete remo	ove eviden	ce card. Attach the ca	ard to the completed checklist.				
					ks are complete or type Q and				
	iter on keyboard.	`			1 31 2				
	100								
SIGNATI	RE								
*WITH (	R WITHOUT VAPO	R RECIRCULATI	ON AND	WITH OR WITHOU	JT KEYBOARD				

DHS/DSLS/Form C116 (Rev. 7-93)

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/14/2004 11:35 COPY NUMBER 1 OF 2 OFFICER BADGE=0000974 N.C.I.C NO.=AZ0071500 BRAC VALUE TIME TEST AIR BLANK .000 11:34 . 100 CAL, CHECK 11:35 .000 11:35 AIR BLANK OPERATOR: ADDITIONAL INFORMATION:

### EXHIBIT Q THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R9-14-404(A)

#### ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

AGENCY GOODYFAR PD QA SPECIALIST WILL
INTOXILYZER SERIAL # (0600/4/6 LOCATION GP) (Print Name)  DATE 1/14 Equ TIME 1/38
Display reads "READY TO START" or "PUSH BUTTON TO START TEST."
DIAGNOSTIC TESTS  DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is 330.  Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".  Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".  Clock time check. Mode selection switch S10 on or keyboard menu selection "E".  Date check. Mode selection switch S11 on or keyboard menu selection "E".
OPERATIONAL TESTS  Alcohol-free subject Test result 0. OCO AC  Error recognition logic system functioning  Invalid test printed  Proper sample recognition system  Invalid sample printed  Deficient sample printed  Completeness of sample purge with collection tube **  Calibration standard 0. OO ACResults: 0. OOO  Instrument operating properly and accurately. YES NO
COMMENTS:
SIGNATURE
O. Alcohol concentration standard collection for subsequent analysis.**
*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD **THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR SAMPLE CAPTURE

DHS/DSLS/Form C117 (Rev. 7-93)

GOODYEAR PD-JUSTICE FACILITY
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
06/14/2004
COPY NUMBER 1 OF 1

11:38

OF.	ļ:	I	CE	Fet	$\mathbb{B}_{!}$	Α	D)	G	Ε.	::::	()	O	()	0	9	7	:: -
Ν.	<u></u> :	ıt	Ι.	ŧ <u></u> :	N	0	n	::::	Α	Z	()	Ö	7	1	5	O	()

DVM TEST

OPERATOR: _		
ADDITIONAL	INFORMATION:	

GOODYEAR PD-JUSTICE FACILITY
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
06/14/2004
COPY NUMBER 1 OF 1

11:38

OF	ļ::··	Ι	[_]	E Fe	BADGE=0	000974	
N.	<u>.</u>	*1	Ţ	, l'''	NO.=AZO	071500	

DISPLAY TEST

OPERATOR:	Andrews Committee of the Committee of th
ADDITIONAL INFORMATION:	***************************************

SN 66-004416 VERSION E 773.23 DATE 06/14/2004 TIME 11:39 COPY NUMBER 1 OF 1

OFFICER BADGE=0000974 N.C.I.C NO.=AZ0071500 ABCDEFGHIJKLMNOPQRSTUVWXYZ0123456789!@#\$ OPERATOR: ADDITIONAL INFORMATION: GOODYEAR PD-JUSTICE FACILITY
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
06/14/2004
COPY NUMBER 1 OF 1

11:39

OF	-	I	CE	F.	В	Α	D	0	E	===	0	0	0	0	9	7	4
Ν.	ı":		Τ.	r":	M	Ü		:::	4	Z	Ō	Ö	7	1	5	Ö	Ó

TIME, DATE, LOC SET

OPERATOR: _		
ADDITIONAL	INFORMATION:	

GOODYEAR PD-JUSTICE FACILITY
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
06/14/2004 11:40
COPY NUMBER 1 OF 2

SUB NAME=TEST, SQAP-AF,
SUB DOB =11/11/11
SUB SEX =MALE WEIGHT =001
OFFICER BADGE=0000974
N.C.I.C NO.=AZ0071500

TEST BRAC VALUE TIME
AIR BLANK .000 11:40
SUBJECT TEST .000 11:40
AIR BLANK .000 11:40
OPERATOR:
ADDITIONAL INFORMATION:

SN 66-004416 DATE 06/14/2004 COPY NUMBER 1 OF 2 VERSION E 773.23 TIME 11:41

INVALID TEST

SUB NAME=TEST,SQAP-PB,
SUB DOB =11/11/11
SUB SEX =MALE WEIGHT =001
OFFICER =WILL,W.C
OFFICER BADGE=0000974
N.C.I.C NO.=AZ0071500

PRESSED START TEST AT IMPROPER TIME

OPE	RA"	FOR	35 14	****					***********		 			
ADD	I T	(ON	AL.	IN	FOF	MA:	TI(	JN:		•••••	 			***************************************
				••••••••							 			
•••••••••••••••••••••••••••••••••••••••	••••••		***********		***************************************	*********	#\$### - T- a T- Y		*********	***************	 ***************************************	***************************************	***************************************	***************************************

SN 66-004416 DATE 06/14/2004 COPY NUMBER 1 OF 2 VERSION E 773.23 TIME 11:42

INVALID TEST

SUB NAME=TEST,SQAP-BWT,
SUB DOB =11/11/11
SUB SEX =MALE WEIGHT =001
OFFICER =WILL,W.C
OFFICER BADGE=0000974
N.C.I.C NO.=AZ0071500

SAMPLE INTRODUCED AT IMPROPER TIME.

OPERATOR:		
ADDITIONAL	INFORMATION:	

/

SN 66-004416 VERSION E 773.23 DATE 06/14/2004 TIME 11:43 COPY NUMBER 1 OF 2

INVALID TEST

SUB NAME=TEST, SQAP-RFI, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =001 OFFICER =WILL, W, C OFFICER BADGE=0000974 N.C.I.C NO.=AZ0071500

INTIEL (ED - EF I
OPERATOR:
ADDITIONAL INFORMATION:

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/14/2004 11:44 COPY NUMBER 1 OF 2 SUB NAME=TEST, SQAP-MOA, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =001 OFFICER =WILL, W, C OFFICER BADGE=0000974 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK .000 11:44 INVALID SAMPLE .MDA 11:44 AIR BLANK .000 11:44 BEGIN NEW DEPRIVATION PERIOD OFERATOR: ADDITIONAL INFORMATION:

GOODYEAR PD-JUSTICE FACILITY INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 06/14/2004 11:49 COPY NUMBER 1 OF 2 SUB NAME=TEST, SQAP-DS, SUB DOB =11/11/11 SUB SEX =MALE WEIGHT =001 OFFICER =WILL, W, C OFFICER BADGE=0000974 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK .000 11:46 \*SUBJECT TEST .097 11:48 AIR BLANK .000 11:49 \* DEFICIENT SAMPLE - VALUE PRINTED WAS HIGHEST OBTAINED. OPERATOR: ADDITIONAL INFORMATION:

### EXHIBIT P THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C R9-14-404(A)

#### ARIZONA DEPARTMENT OF HEALTH SERVICES

### STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

STANDARD CALIBRATION CHECK PROCEDURE

Agency	′	(500DY=AR 8D) Date 4(Ko/01/ Time /15						
Intoxily	zer S	Serial # 106-0044b Location 69						
QA Spe	ciali	st Will						
(Print Name)								
$\langle \mathcal{W} \rangle$	1.	Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the						
i17.		simulator. Insure that a tight seal has been made. Standard value: 0. AC						
(W)	2.	•						
$(\mathcal{U})$	3.	Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or						
. /	selecting "C" on keyboard menu.							
$(\mathcal{V})$	4.	Attach simulator to the simulator entrance port on the Intoxilyzer.						
$(\mathcal{W})$	5.	Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".						
(V&)	6.	Push Start Test button or press enter on keyboard.						
$(\mathcal{W})$	7.	Insert card in response to display.						
$(\mathcal{W})$	8.	Air blank completed.						
(W)	9.	Calibration check completed. Test results 0.099AC						
(W)	10.	Air blank completed.						
$(\mathcal{W})$	11.	When display reads Test Complete remove evidence card. Attach the card to the completed checklist.						
$(\mathcal{W})$	12.	Return mode selection switch #9 to off position after all calibration checks are complete or type Q and						
		er on keyboard.						
		A						
SIGNA	TUR	E						
*WITH	*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD							

DHS/DSLS/Form C116 (Rev. 7-93)

GOODYEAR PD-JUSTICE FACILITY
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
06/14/2004 11:51
COPY NUMBER 1 OF 2

OFFICER BADGE=0000974
N.C.I.C NO.=AZ0071500

TEST BRAC VALUE TIME
AIR BLANK .000 11:50
CAL. CHECK .099 11:51
AIR BLANK .000 11:51

OPERATOR:

ADDITIONAL INFORMATION: