# EXHIBIT P THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C R9-14-404(A)

#### ARIZONA DEPARTMENT OF HEALTH SERVICES

## STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

STANDARD CALIBRATION CHECK PROCEDURE

Agency		2000 YEAR +0 Date 9-9-02 Time 1635					
Intoxilyzer Serial # (66 604416 Location 600							
QA Specialist WW							
(Print Name)							
<b>(W</b> )	1.	Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0. / CCC					
$(\mathcal{M})$	2.	Turn on the simulator and allow the temperature to reach 34°C±.2°C.					
( <b>1</b> )	3.	Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.					
$(\mathbf{w})$	4.	Attach simulator to the simulator entrance port on the Intoxilyzer.					
(W)	5.	Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".					
( <b>W</b> )	6.	Push Start Test button or press enter on keyboard.					
( <b>W</b> )	7.	Insert card in response to display.					
( <b>W</b> )	8.	Air blank completed.					
( <b>W</b> )	9.	Calibration check completed. Test results OPI AC					
( <b>W</b> )		Air blank completed.					
$(\mathbf{w})$		When display reads Test Complete remove evidence card. Attach the card to the completed checklist.					
(W		Return mode selection switch #9 to off position after all calibration checks are complete or type Q and					
	ente	er on keyboard					
SIGNATURE							
*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD							

DHS/DSLS/Form C116 (Rev. 7-93)

GOODYEAR POLICE DEPT.
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004416
09/09/2002 16:36
COPY NUMBER 1 OF 2

OFFICER BADGE=0000974
N.C.I.C NO.=AZ0071500

TEST BRAC VALUE TIME
AIR BLANK .000 16:35
CAL. CHECK .093 16:36
AIR BLANK .000 16:36

OPERATOR:

ADDITIONAL INFORMATION:

## EXHIBIT P THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C R9-14-404(A)

#### ARIZONA DEPARTMENT OF HEALTH SERVICES

## STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

STANDARD CALIBRATION CHECK PROCEDURE

Agency		bodyear PD Date 9-11-02 Time 0250				
	$\overline{}$	erial # 66 - 00 4416 Location 6PD				
		st Will				
(Print Name)						
( N	1.	Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0.7 to AC				
(1)	2.	Turn on the simulator and allow the temperature to reach 34°C ±.2°C.				
( <b>%</b> )	3.	Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.				
( <b>~</b> )		Attach simulator to the simulator entrance port on the Intoxilyzer.				
(1)	5.	Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".				
$(\mathcal{Y})$	6.	Push Start Test button or press enter on keyboard.				
$(\mathcal{N})$	7.	Insert card in response to display.				
$\mathcal{W}$		Air blank completed.				
$(\mathcal{N})$	9.	Calibration check completed. Test results 0.893AC				
$\langle \hat{\mathcal{W}} \rangle$	10.	Air blank completed.				
(W)		When display reads Test Complete remove evidence card. Attach the card to the completed checklist.				
(VV)	12.	Return mode selection switch #9 to off position after all calibration checks are complete or type Q and				
•	ente	r on keyboard				
SIGNA	TUR	E				
*WITH	OR	WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD				
** * * * * * * * * * * * * * * * * * * *		THE THE STATE OF T				
DHS/DS	SLS/I	Form C116 (Rev. 7-93)				

GOODYEAR POLICINTOXILYZER - AZ MODEL 5000 09/11/2002 COPY NUMBER 1	07:51					
OFFICER BADGE=0000974 N.C.I.C NO.=AZ0071500						
AIR BLANK CAL. CHECK AIR BLANK	RAC VALUE .000 .093 .000	TIME 07:50 07:50 07:51				
ADDITIONAL INFORMATION:						
		A				

## EXHIBIT P THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C R9-14-404(A)

#### ARIZONA DEPARTMENT OF HEALTH SERVICES

### STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 5000\*

STANDARD CALIBRATION CHECK PROCEDURE

Agency	ے) ر	Deo Dye AR +D Date 9-11-02 Time 0 75 9				
Intoxily	zer S	Serial # 66 - 0416 Location 690				
QA Specialist WILL						
(Print Name)						
$(\mathbf{W})$	l.	Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0.				
( <b>W</b> )	2	Turn on the simulator and allow the temperature to reach 34°C±.2°C.				
(V)	3.					
(ك	4.	Attach simulator to the simulator entrance port on the Intoxilyzer.				
( <b>W</b> )	5.	Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".				
$(\mathcal{U})$	6.	Push Start Test button or press enter on keyboard.				
$(\mathcal{V})$	7.	Insert card in response to display.				
$(\mathcal{N})$	8.	Air blank completed.				
$\langle \mathcal{N} \rangle$	9.	Calibration check completed. Test results 0.90 AC				
$(\mathcal{N})$	10.	Air blank completed.				
	11.	When display reads Test Complete remove evidence card. Attach the card to the completed checklist.				
	12.	Return mode selection switch #9 to off position after all calibration checks are complete or type Q and				
• •	ente	er on keyboard				
SIGNA						
*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD						

DHS/DSLS/Form C116 (Rev. 7-93)

INTOXILYZER - ALCOHOL ANALYZER AZ MODEL 5000 SN 66-004416 09/11/2002 08:00 COPY NUMBER 1 OF 2 OFFICER BADGE=0000974 N.C.I.C NO.=AZ0071500 TEST BRAC VALUE TIME AIR BLANK .000 07:59 .097 CAL. CHECK 08:00 AIR BLANK .000 08:00 OPERATOR: ADDITIONAL INFORMATION: SOCUTION CHANGE ACS 700701B

GOODYEAR POLICE DEPT.