### EXHIBIT G-2 THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

#### ARIZONA DEPARTMENT OF PUBLIC SAFETY

### STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE OA SPECIALIST DATE 12-20-06 LOCATION GOODYEAR INTOXILYZER SERIAL# 80-000312 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution AC. Pour a standard alcohol concentration solution \_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
 Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER. 4. Air blank completed. 5. Calibration check completed. Test results 0. 10 6. Air blank completed. -7. Remove printed record, Attach the record to the completed checklist. **SIGNATURE** 

DPS Form Exh G-2 (Rev 05+01)

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

12/20/2006

15:13:30

Standard Lot#: 618601I

Last Changed By: J. ROGERS #1123

QAS: J. ROGERS #1123

GOODYEAR PD

-----

Test	g/210L	Time
Air Blank	0.000	15:13:57
Cal Check	(0.101)	15:14:16
Air Blank	0.000	15:14:46

#112

Operator s Signature

## EXHIBIT G-4 THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

#### ARIZONA DEPARTMENT OF PUBLIC SAFETY

## STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE Display Reads "PUSH BUTTON TO START" DIAGNOSTIC TESTS

Al. Clock time check. 2. Date check. OPERATIONAL TESTS 1. Alcohol-free subject test result 0. O AC. 2. Error recognition logic system functioning. Not a Successfully Completed Test Sequence printed -3. Proper sample recognition system. Not a Successfully Completed Test Sequence printed Deficient sample printed. 4. Standard Calibration Check standard 0. 100 AC. Result 0. 100 AC. Instrument is operating properly and accurately. YES **COMMENTS** -#1123 SIGNATURE\_ DPS Form Exh G-4 (Rev 05/01)

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

12/20/2006 15:15:58

Standard Lot#: 618601I Last Changed By: J. ROGERS	#1123
QAS: J. ROGERS GOODYEAR PD	#1123
Operator: J. ROGERS GOODYEAR PD	#1123

Subject: TEST, TEST DOB: 01/01/0101

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	15:16:49
Diagnostic Test	Pass	15:17:17
Air Blank	0.000	15:17:44
0.100 Cal Check	0.101	15:18:03
Air Blank	0.000	15:18:32
Subject Test (	(0.000)	15:18:51
Air Blank	0.000	15:19:18
Five Minute Wait(	(ABT* )	
Air Blank	0.000	15:19:56

<sup>\*</sup>Sequence Aborted

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

12/20/2006 15:20:22

Standard Lot#: 618601I Last Changed By: J. ROGERS	#1123
QAS: J. ROGERS GOODYEAR PD	#1123
Operator: J. ROGERS GOODYEAR PD	#1123

Subject: TEST, TEST DOB: 01/01/0101

Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	(IPS*)	15:20:58
Air Blank	0.000	15:21:27

<sup>\*</sup>Improper Sample

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

12/20/2006 15:21:48

Standard Lot#: 618601I

Last Changed By: J. ROGERS #1123

QAS: J. ROGERS #1123

GOODYEAR PD

Operator: J. ROGERS #1123

GOODYEAR PD

Subject: TEST, TEST DOB: 01/01/0101

Sex: M Weight: 100

15 Minute Deprivation Period? Yes

 Test
 g/210L
 Time

 Air Blank
 RFI\*
 15:22:28

 Air Blank
 0.000
 15:22:57

<sup>\*</sup>RFI Detect

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

12/20/2006 15:23:22

Standard Lot#: 618601I Last Changed By: J. ROGERS	#1123
QAS: J. ROGERS GOODYEAR PD	#1123
Operator: J. ROGERS GOODYEAR PD	#1123

Subject: TEST, TEST DOB: 01/01/0101

Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	15:24:10
Diagnostic Test	Pass	15:24:38
Air Blank	0.000	15:25:05
0.100 Cal Check	0.101	15:25:24
Air Blank	0.000	15:25:54
Subject Test	0.015*	15:29:02
Air Blank	0.000	15:29:31
Air Blank	0.000	15:29:58
Subject Test (	XXX**	15:30:16
Air Blank	0.000	15:30:50

<sup>\*</sup>Deficient Sample

Begin new deprivation period

<sup>\*\*</sup>Invalid Sample

# EXHIBIT G-2 THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

#### ARIZONA DEPARTMENT OF PUBLIC SAFETY

# STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

#### STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST T. ROGERS	AGENCY GOODYEAR P.D.
QA SPECIALIST	AGENCY GOOTETT
DATE 12-20-06	TIME 1>31
INTOXILYZER SERIAL# 80-000312	LOCATION GOODYEAR P.D.
1. Ensure that gas tank is attached to instrument and	contains a standard alcohol concentration solution ACOR
Pour a standard alcohol concentration solution simulator. Ensure that a tight seal has been made. 0.2°C  2. Intoxilyzer 8000 display reads "PUSH BUTTON"	AC, into a clean dry simulator and assemble the Turn on the simulator and allow temperature to reach 34°C ± TO START" dry control test or "W" for wet control test. After selection isAC.
√ / /	

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

12/20/2006

15:31:09

-----

Standard Lot#: 618601I

Last Changed By: J. ROGERS #1123

QAS: J. ROGERS #1123

GOODYEAR PD

------

Test	:	g/210L	Time
Air	Blank	0.000	15:31:37
Cal	Check	0.100	15:31:55
Air	Blank	0.000	15:32:25

Operator's Signature